

TECHNICAL DATA SHEET**Catalog # MS507****ATP synthase subunit alpha monoclonal antibody**

Rev.3

LOT #:**COMPONENTS:** 100 µg monoclonal antibody**APPLICATIONS:** Western blotting, Immunocytochemistry, Flow Cytometry**CLONE ID OF MONOCLONAL ANTIBODY (mAb):** 15H4C4**SPECIES CROSS-REACTIVITY:** human, bovine, mouse, monkey, rat, *Drosophila melanogaster*, *C. elegans*, and oyster**HOST SPECIES AND ISOTYPE:** Mouse IgG2b, k**IMMUNOGEN:** Bovine Complex V**CONCENTRATION:** 1 mg/mL in HEPES-Buffered Saline (HBS) with 0.02% azide as a preservative.**SUGGESTED WORKING CONCENTRATION:**
1 µg/mL for Western blotting
1-2 µg/mL for Immunocytochemistry
1 µg/mL for Flow Cytometry**mAb PURITY:** Near homogeneity as judged by SDS-PAGE. The antibody was produced *in vitro* using hybridomas grown in serum-free medium, and then purified by biochemical fractionation.**STORAGE CONDITIONS:** Store at 4°C. Do not freeze.**COUNTRY OF ORIGIN:** USA**BACKGROUND:**

Complex V, also called F_1F_0 ATPase or ATP synthase, is responsible for ATP production in oxidative phosphorylation and can work in reverse as a proton pumping ATPase. The enzyme was thought to be localized exclusively to mitochondria. However, it has recently been identified on the plasma membrane of several cell types including hepatocytes where it functions as the HDL receptor, on endothelial cells where it may act as the angiotensin receptor, and on the surface of cancer cells.

The enzyme in mammals is composed of 17 subunits, five of which make up the easily detached F_1 . The remainder subunits are components of two stalk domains and the proton pumping F_0 part of the machinery. Two of the subunits of the F_0 part are encoded on mitochondrial DNA while the other subunits are nuclear encoded. Mutations in the mitochondrial-encoded subunits of ATP synthase (Complex V) cause OXPHOS disease.

Note: This product is for research purposes only. It is not to be used in humans or for diagnostic purposes.

WARRANTY – MitoSciences Inc warrants that products will perform as indicated in the published Product Information Sheet and Product Technical Specifications for 6 months from date of purchase when stored according to specifications and when used consistent with recommended protocols. If you experience results which materially differ from those described, send evidence of the non-performing product for replacement of the original product purchased or a credit toward any other of the company's products or services.